Permanent Crown

EN	Instructions for Use	NL	Gebruiksaanwijzing
DE	Gebrauchsanweisung	BG	Инструкция за употреба
FR	Notice d'utilisation	SV	Bruksanvisning
ES	Instrucciones de uso	PL	Instrukcja używania
IT	Istruzioni per l'uso	SL	Navodila za uporabo
PT	Instruções de utilização		

10/ 13/ 2020 PRNT-106964-02

Permanent Crown

Resin for 3D printing of permanent single crowns, inlays, onlays and veneers.

1. Intended Use / Indication

Permanent Crown is a light-curing, free-flowing plastic based on methacrylic acid esters for the production of permanent single crowns, inlays, onlays and veneers.

2. Contraindications

Known allergy to one or more ingredients. In cases of doubt, the allergy should be clarified and ruled out based on a specific test prior to the application of this product. Permanent Crown should not be used for purposes other than the production of permanent single crowns, inlays, onlays and veneers. Any deviation from these instructions for use can have negative effects on the chemical and physical quality of plastics made from Permanent Crown.

3. Safety instructions

Permanent Crown is produced and tested according to the most stringent quality standards. In order to ensure optimum further processing, please read the information contained in the instructions for use carefully. The improper use of Permanent Crown and failure to follow information can have a detrimental effect on the quality of objects produced from Permanent Crown. Nitrile gloves, safety goggles and a coat must be worn as a means of protection when handling the resin and the plastic that has not been post-cured. Conventional medical gloves do not offer any lasting protection against the sensitizing effect of methacrylates. If the product comes into contact with the glove, take the glove off and discard it, wash your hands immediately with water and soap and put on a new glove. Consult a physician in the event of an allergic reaction.

The safety and care instructions set down in the Permanent Crown instructions for use and safety data sheet shall apply to the handling of liquid resin and printed objects that have not been post-cured (objects in the “green condition”). A dust mask must be worn too due to potential dust formation while the printed objects are being processed.

It is prohibited to use plastic parts made of Permanent Crown as auxiliary equipment for food and drinks applications.

4. Side effects and precautions

PRECAUTIONS / PROTECTION

It is essential that protective clothing be worn when handling Permanent Crown. Safety goggles and nitrile gloves must be used. However, we cannot completely rule out the possibility of personal reactions to individual components in isolated cases. In such cases, the respective user should discontinue use of Permanent Crown. If intolerances or allergic reactions occur when it comes into contact with the patient, discontinue use of the material.

WARNING

Contains:

Esterification products of 4,4'-isopropylidiphenol, ethoxylated and 2-methylprop-2enoic acid. Silanized dental glass, methyl benzoylformate, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide. Total content of inorganic fillers (particle size 0.7 µm) is 30–50% by mass.

Hazard statements as per MSDS

- Causes skin irritation.
- May cause an allergic skin reaction.
- Causes serious eye irritation.
- May cause respiratory irritation.
- May cause long harmful effects to aquatic life.

Precautionary statements as per MSDS

- Avoid breathing mist / vapours / spray.
- Avoid release to the environment.
- Wear protective gloves / protective clothing / eye protection / face protection.
- Call a POISON CENTER / doctor if you feel unwell.
- If skin irritation or rash occurs: Get medical advice / attention.
- If eye irritation persists: Get medical advice / attention.
- Dispose of contents / container as per local and national regulations.

5. General information on handling

Permanent Crown is supplied in colours according to the VITA* classical shades.

STORAGE

Permanent Crown must be stored in the original cartridge at room temperature approx. 22 °C in a dark, dry place. It must be ensured that the temperature does not drop below +4 °C and does not exceed +28 °C! The minimum shelf life date printed on the product must be observed. Do not use resin cartridges past the

*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

expiration date on the cartridge label. Perfect processing cannot be guaranteed if materials which have exceeded their minimum shelf life date are used. The completely cured print objects must be stored at room temperature and protected from sources of bright light.

6. Processing

Permanent Crown is one of the system components in the Formlabs 3D print system and has been optimized for use in the SLA 3D printer – Form 3B.

NOTICE:

SINGLE CROWNS, INLAYS, ONLAYS AND VENEERS

Minimum wall thicknesses anterior teeth	1.0 mm
---	--------

Minimal wall thicknesses posterior teeth	1.0 mm
--	--------

Minimal wall thicknesses cervical	1.0 mm
-----------------------------------	--------

All accessories used must be designated for Permanent Crown only. For full compliance and biocompatibility of restorations from Permanent Crown, the process requires a dedicated resin tank, Stainless Steel Build Platform, Form Wash, and finishing kit which should not be mixed with any other resins.

The ideal working temperature of Permanent Crown is in the temperature range between 18 °C and 28 °C . **Before using a new cartridge of Permanent Crown resin, the material has to be shaken well!** For further processing – selecting the resin, setting up the print job – as part of the printing process, follow the Formlabs Form 3B printer instructions for use. Permanent Crown in the resin tank does not need to be mixed prior to printing, up to a month of non-use, as solids (fillers) in the material do not settle.

PRINTING:

- a. Insert resin cartridge, Stainless Steel Build Platform, and the compatible resin tank into the Formlabs 3D printer.
- b. Prepare a print job using PreForm software. Import a dental restoration STL file. Orient files horizontally with the occlusal plane facing the build platform and generate supports. For recommendations on print orientation and support placement, see the detailed application guide at support.formlabs.com.
- c. Send print job to printer. Begin print by selecting print job from print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete the print.

Subsequent processing

On completion of printing, the print objects are released from the build platform using the spatula supplied. The print objects should be cleaned using the Form Wash as follows.

WASHING:

- a. Remove the build platform with the printed parts from the printer and place it into a Form Wash with clean IPA ($\geq 99\%$).
- b. Set the Form Wash for 3 minutes. Make sure the parts are fully submerged in IPA when washing. Exceeding wash duration may affect dimensional accuracy and performance of printed parts over time.
- c. Once parts are washed, inspect parts for any remaining uncured resin on print surfaces. Use an IPA-filled spray bottle to remove any remaining resin on the printed parts and in between supports and rafts. An IPA-soaked brush may be used to help with uncured resin removal.

DRYING: Use compressed air to dry the parts. The dried parts will have a white, powdery coating on the print surfaces. If any wet, uncured resin is still present after drying, reuse the IPA-soaked brush and IPA-filled squeeze bottle to remove uncured resin and dry again.

REMOVAL: After drying, remove printed parts from the build platform by wedging the part removal tool under the print raft and rotating the tool. For detailed techniques, see the support material on dental.formlabs.com.

POST CURE: The final properties of the print object depend on the post-curing process. These are achieved by light polymerization with the Formlabs Form Cure. To maintain dimensional accuracy and biocompatibility, post curing must be done in two steps:

1. With the raft and supports still intact, place printed parts in the Form Cure with the raft side down. Form Cure parts at 60°C for 20 minutes.
 - a. Remove parts from Form Cure and use a handpiece with a cutting disc or side cutters to separate supports and raft from printed parts.
 - b. Carefully sandblast the printed part surfaces to remove the white coating using a sandblaster with glass bead blasting material 50 μ m, (e.g. Perlablast® micro*, BEGO, REF 46092/54302) at a maximum blasting pressure of 1.5 bar.
 - c. Check printed part fit on the model and finish and contour as needed using a dental handpiece and carbide rotary burs or diamond grinding stones.
2. Place finished parts in the Form Cure and post cure again at 60°C for an additional 20 minutes, with the occlusal plane facing upwards. For best results

*This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

and full compliance, printed parts must be post cured twice as recommended. Any deviations from this protocol, may lead to unsatisfactory fit and affect biocompatibility.

Notice: The times given only apply to regularly maintained equipment that guarantees a corresponding light intensity. Inspect the dental restorations for cracks. Discard if any cracks or damage are detected.

INDIVIDUALIZATION: Optional step: Individualization of the post-cured objects is possible with composite stains and is the responsibility of the user. It may affect the color result. Follow manufacturer's instructions. If this step is not considered, follow the next point.

POLISHING: Polish the surface of the objects with pumice stone and polishing compound

- a. Avoid overheating the parts during polishing.
- b. Optimal surface quality is achieved by polishing after post-curing.
- c. If rough marks are left, sand and polish to smooth surfaces to improve patient comfort.

7. Cleaning in the dental laboratory and dental practice

Fully cured objects made from Permanent Crown can be easily cleaned and disinfected. Steam cleaning (e. g., with Triton SLA*) is possible. Disinfection in the immersion bath (e. g. ethanol 96 % or MD 520* impression disinfectant, Dürr Dental Co.) is also possible. Follow manufacturer's instructions.

8. Notice for practitioners

Restorations can undergo high-gloss polishing with composite polishers commonly used in dental practice. The finished restorations can be attached using self-adhesive cements (e. g. RelyX Unicem*, 3M Espe*) or composite cement with a primer (e. g. Variolink Esthetic DC* and Monobond Plus*, Ivoclar Vivadent*). Observe the instructions for use of the luting agent.

9. Disposal

The cured, separated material (base plate, support structure) can no longer be used. Cured material can be disposed of as domestic waste. Unused resin, contaminated packaging or IPA used for cleaning with resin residues must be disposed of via the local waste disposal authority or a hazardous waste collection point stating the safety data sheet.

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Herbst GmbH & Co. KG**