GQ QUATTRODISC NEM SOFT

GEBRAUCHSINFORMATION INSTRUCTIONS FOR USE

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EN

DESCRIPTION

Beryllium- and Cadmium-free cobalt base blanks are indicated as Type 4 for the manufacture of fixed and removable dentures for crowns and bridge frameworks as well as bridge and implant superstructures according to DIN EN ISO 22674 and EN ISO 9693. The processing is carried out using CAD/CAM technology with suitable milling machines (Please note manufacturer's specifications.).

SCOPE OF APPLICATION

- Single crowns
- Bridges in every span width for veneering with dental ceramics
- Cone and telescope technology
- Implant construction
- Primary and secondary parts
- Full cast crowns and full cast bridges

MATERIAL CHARACTERISTICS

Type of alloy	4
Density	8.8 g/cm ³
Vickers hardness	285 HV 10/30
CTE value (25–500 °C)	14.4 × 10 ⁻⁶ K ⁻¹
Oxidation without vacuum	1,000 °C
Mass content of the components (in % by weight)*	Co 63, Cr 29, Mo 6, Mn x, Nb x, Si x, Fe x
Tensile strength	680 MPa
0.2% proof stress	400 MPa
Elongation at break	9 %
Young's modulus	250 GPa

* x = <1 %

PROCESSING INFORMATION

Fabrication

Insert the disc into the machine according to the operating instructions of the CAM system and process the disc using tungsten carbide burrs suitable for CoCr. When finished, remove the blank with the frameworks produced

Separating the frameworks

The milled restorations can be cut out of the blank using cross-cut tungsten carbide burrs or separating discs suitable for CoCr alloys.

Finishing/Cleaning

Machine the separated frameworks with sharp tungsten carbide burrs without high pressure and in one direction only. Avoid sharp edges and interfacing areas. For the desired veneering, the surface must be blasted with aluminium oxide (110–180 µm) at a pressure of 3–4 bar, thoroughly evaporated and then degreased with ethanol.

ATTENTION: Metal dust is harmful to health. During finishing and stripping, use suitable suction and a respirator mask of type FFP3-EN-2021

Oxidising

An oxide firing is not necessary with CoCrMo₆. If, however, an oxide firing (rapid heating to 1.000 °C without vacuum, cooling down immediately after reaching the temperature) is carried out for visual inspection of the framework conditioning, the oxide must subsequently be stripped off with single-use abrasive aluminium oxide (110 –180 µm) and subsequently cleaned.

Veneering with ceramics

All conventional ceramic materials with an appropriate CTE value are suitable for ceramic veneering. The processing instructions recommended by the ceramic manufacturer must be observed, particularly with regard to the cooling conditions after ceramic firing.

Joining (soldering/lasering)

For soldering, a solder and flux suitable for CoCr alloys is recommended. Dimension the soldering block as small as possible and heat it initially at about 600 °C for 10 minutes in the oven. Soldering gaps should be between 0.05 and 0.2 mm to avoid solder flowing through solder joints. For laserwelding prior firing, a special laserwelding wire is recommended.

Polishing

After ceramic firing or soldering, existing oxides/flux residues must be removed and the framework finished to a high gloss using suitable grinding/polishing instruments and polishing pastes.

NOTES

Contraindications

- In case of known hypersensitivity (allergy) to constituents of the alloy.
- All applications not listed under Indications are contraindicated.

Side effects

In individual cases, hypersensitivity reactions to constituents of the alloy used cannot be ruled out. The product should not be used in case of known allergies to the constituents of the alloy.

Interactions

On contact with or in close proximity to alloy of a different type, galvanic effects can occur. In event of sustained local discomfort, other materials must be used for the denture.

SAFETY INSTRUCTIONS

- When designing the prosthetic dentures, ensure adequate wall thickness and connector cross sections.
- These medical devices may only be used and applied by qualified staff in a dental laboratory.
- Please also observe the material safety data sheet.
- Metal dust is harmful to health. During finishing and stripping, use suitable suction and a respirator mask of type FFP3 EN-2021
- Ensure batch traceability even after the milling process

DISPOSAL

Dispose of contents in accordance with local, regional, national or international regulations.

Gold Quadrat GmbH

Büttnerstr. 13
30165 Hannover
Germany

T +49 (0)511 449897-0
info@goldquadrat.de

www.goldquadrat.de

Wieland Edelmetalle GmbH

Schwenninger Str. 13
75179 Pforzheim
Germany

T +49 (0)7231 1393-0

www.wieland-edelmetalle.de

CE
0483

GQ2111/02.24