GEBRAUCHSINFORMATION INSTRUCTIONS FOR USE

GQ METALLICAST NP

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EN

DESCRIPTION

GQ MetalliCast NP is a cobalt-based alloy according to EN ISO 22674 for restorations with thin cross-sections that are subject to high forces, e.g. thin veneered single crowns, fixed multi-unit bridges or bridges with small cross-sections, bars, implant-supported superstructures and removable dentures.

GQ MetalliCast NP is characterised by its high flowability and very good melting behaviour and allows the thinnest castings up to 0.3 mm.

GQ MetalliCast NP has a smooth and compact surface. A reduced oxide formation, even during repeated and long firing processes in the ceramic furnace, is given with this alloy. Thanks to its ideal thermal expansion coefficient, this alloy is suitable for all types of last generation ceramics.

The alloy is highly resistant to corrosion and does not contain nickel.

SCOPE OF APPLICATION

- Single crowns
- Bridges in every span width for veneering with dental ceramics
- Cone and telescope technology
- Primary and secondary parts
- Implant construction
- Full cast crowns
- Full cast bridges

MATERIAL CHARACTERISTICS

Alloy type	4
Color	white
Density	8.4 g/cm ³
Solidus-liquidus temperature	1,253–1,304 °C
Melting point	1,360 °C
Vickers hardness	324 HV10
CTE value (25–500 °C)	14.1 × 10 ⁻⁶ K ⁻¹
CTE value (25–600 °C)	14.5 × 10 ⁻⁶ K ⁻¹
Mass content of the components (in % by weight)*	Co 63.0; Cr 28.0; Nb 4.0; W 3.0; Si 1.0; Mn x; Fe x
0.2% proof stress	475 MPa
Elongation at break	3 %
Young's modulus	194 GPa
Ions release after 7 days	0.8 µg/cm ²
Maximum baking temperature	980 °C

* x = <1 %

PROCESSING INFORMATION

Modelling

Model with a minimum thickness of 0.3 mm for simple structures and 0.5 mm in very extensive cases or in patients with bruxism. Avoid acute angles. Ensure a reduced anatomical shape for copings intended for ceramic veneering.

Casting channels

Sprue channels for direct sprueing: Ø 3 mm, feed channels Ø ≥ 6 mm
Sprueing channels for indirect sprueing with stabilising bar: Ø 3 mm,
Beam Ø ≥ 5 mm, feed channels Ø ≥ 6 mm

Casting

Pour GQ MatalliCast NP in a clean and preheated ceramic crucible.

Each alloy must have its own specific crucible, do not preheat the alloy and do not use flux.

A preheating temperature of 850–950°C must be achieved. A holding time depends on the specifications of the investment material manufacturer. Ensure that the investment ring is heated evenly and completely.

Induction casting

Start casting when the alloy is close to breaking up the oxide layer.

Casting with burner

Place the alloy in the preheated ceramic crucible and heat it up with circular movements. When the alloy is melted, operate the centrifuge. Only burners with multiple flames can be used. Do not use fluxes.

Guide values for flame regulation:

- acetylene 0.4 bar, oxygen 2 bar
- propane 0.2 bar, oxygen 2 bar
- methane pressure, oxygen 2 bar

PROCESSING

Allow the casting to cool at ambient temperature, remove the investment material and sandblast with aluminium oxide from 110–150 µm at a pressure of 3–4 bar.

Homogenisation and degassing of the cast workpiece at 1,000°C for 10 minutes in the atmosphere is recommended.

For further machining of the cast workpieces, milling cutters made of tungsten carbide are recommended.

After machining the piece, sandblast it with alumina from 110–150 µm at a pressure of max. 3–4 bar and finally clean it with steam. After cleaning, the workpiece should not be touched and should be held only with clamping pliers.

OXIDATION

Oxidise from 950°C to 980°C for 5 minutes. After firing, sandblast the oxide layer again with aluminium oxide from 110–150µm at a pressure of 2.5–4 bar.

Then evaporate and check that the surface is uniformly grey.

REUSE OF SPRUES

The best results are obtained with pure GQ MetalliCast NP. However, it is possible to reuse the cast material once, provided that the same amount of the new material is added and that it comes from the same batch.

CERAMIC LAYERING

You can use a bonding. Proceed with the application of the ceramic according to the manufacturer's instructions.

Long-term cooling is recommended.

All ceramic materials suitable for cobalt-based alloys, such as Noritake EX-3 or VITA VM 13, can be used for ceramic veneering.

SOLDERING

For soldering we recommend the use of our alloy GQ Solder Stick.

5 SAFETY INSTRUCTIONS

- These medical products may only be used and applied in the dental laboratory by qualified personnel!
- In rare cases, the cobalt-chrome alloy can cause skin inflammation when touched by the susceptible persons. A contact test (patch test) is recommended.
- The dust and smoke from the metals can be harmful to health. During melting and polishing, use an appropriate exhaust system.
- Before using the denture, check if there are other metallic dental restorations in the patient's oral cavity. Galvanic effects may occur if different metals are present.
- It should be noted that the material used may have an influence on radiological examinations, as it is radioopaque.

The product is not sterile.

Storage temperature -30°C to 50°C

CAUTION: Metal vapours and metal dusts are hazardous to health when inhaled; therefore always use a suction device and/or suitable protective mask!

6 DISPOSAL

The machining residues must be disposed of as special waste in accordance with the guidelines 2008/98/EC on waste and 94/62/EC on packaging and packaging waste in accordance with the relevant national rules be disposed of.

Gold Quadrat GmbH

Büttnerstraße 13
30165 Hannover

T +49 (0) 511 449897-0
F +49 (0) 511 449897-44

info@goldquadrat.de

www.goldquadrat.de

CE 0197

GQ1883/06.22