

# KOMBINATIONSTABELLE

VITA 3D-MASTER®	NORITAKE VALUE SHADE		NORITAKE CZR-MASSEN		Enamel	Translucent/Luster <sup>2</sup>
	Body (Dentin)	Opacious Body	Shade Base Stain	Margin		
0M1	NW 0 B	–	SSNW 0	MNW 0	Silky E2	T1/LT1
0M2	NW 0,5 B	–	SSNW 0	MNW 0,5	Silky E2	T1/LT1
0M3	NW 0,5 B	–	SSNW 0,5	MNW 0,5	Silky E2	T1/LT1
1M1	1110 B	OB 1110	SSA1	MA1	E1	T1/LT1
1M2	1120 B	OB 1120	SSA1	MA1	E1	T1/LT1
2L1,5	2015 B	OB 2015	SSB2	MB 2	E1	T1/LT1
2L2,5	2025 B	OB 2025	SSB2	MB2 + MD4 <sup>1</sup>	E1	T1/LT1
2M1	2110 B	OB 2110	SSA1	MA1	E1	T1/LT1
2M2	2120 B	OB 2120	SSA2	MA2	E1	T1/LT1
2M3	2130 B	OB 2130	SSB3	MB3	E1	T1/LT1
2R1,5	2215 B	OB 2215	SSA1	MA2	E1	T1/LT1
2R2,5	2225 B	OB 2225	SSA2	MA2	E1	T1/LT1
3L1,5	3015 B	OB 3015	SSC2	MC2	E2	T1/LT1
3L2,5	3025 B	OB 3025	SSB3	MB3 + MD4 <sup>1</sup>	E2	T1/LT1
3M1	3110 B	OB 3110	SSC1	MC2	E2	T1/LT1
3M2	3120 B	OB 3120	SSA3	MA3	E2	T1/LT1
3M3	3130 B	OB 3130	SSA3,5	MA3,5	E2	T1/LT1
3R1,5	3215 B	OB 3215	SSA2	MA3 + MC2 <sup>1</sup>	E2	T1/LT1
3R2,5	3225 B	OB 3225	SSA3	MA3 + MC2 <sup>1</sup>	E2	T1/LT1
4L1,5	4015 B	OB 4015	SSC2	MA4 + MC2 <sup>1</sup>	E2	T1/LT1
4L2,5	4025 B	OB 4025	SSB3	MA4 + MB3 <sup>1</sup>	E2	T1/LT1
4M1	4110 B	OB 4110	SSC2	MC 2	E2	T1/LT1
4M2	4120 B	OB 4120	SSA4	MA4	E2	T1/LT1
4M3	4130 B	OB 4130	SSA4	MA4	E2	T1/LT1
4R1,5	4215 B	OB 4215	SSA2	MD3 + MC2 <sup>1</sup>	E2	T1/LT1
4R2,5	4225 B	OB 4225	SSA4	MA4	E2	T1/LT1
5M1	5110 B	OB 5110	SSC2	MA4	E2	T1/LT1
5M2	5120 B	OB 5120	SSA4	MA4	E2	T1/LT1
5M3	5130 B	OB 5130	SSA4	MA4	E2	T1/LT1

<sup>1</sup> Mischverhältnis 1:1

<sup>2</sup> OPTION: Es können alternative CZR Translucent/Luster-Massen verwendet werden.